




sinus-Reduction Stent

Self-expandable nitinol stent system




The repositionable
nitinol stent to reduce
the blood flow



REDUCING THE BLOOD FLOW

sinus- Reduction Stent





closed-cell design



radiopaque markers




anti-jump technique



atraumatic soft-tip



repositionable stent




- Zur Erhöhung des Druckgradienten nach TIPS-Anlage
- Repositionierbarer selbstexpandierender Nitinol Stent
- Nutzbare Länge des Einführbestecks: 70 cm



- Augmentation du gradient de pression après mise en place d'un TIPS
- Stent autoexpansible en Nitinol repositionnable
- Longueur système de pose: 70 cm



- To increase the pressure gradient after TIPS procedure
- Repositionable self-expanding nitinol stent
- Length of application device: 70 cm



Bestell-Nr. - Länge des Einführbestecks référence - longueur système de pose order code - length of application device	Aussen extérieure exterior	Taille de l'étranglement waistline	Länge longueur length
10F / 70 cm	Ø / mm	Ø / mm	mm
7542-6540*	12	4	40
7544-6540	14	4	40
7552-6540*	12	5	40
7554-6540	14	5	40
7564-6540*	14	6	40

* Mindestbestellmenge / Lieferzeit auf Anfrage
 * quantité minimum de commande / délai de livraison sur demande
 * minimum order quantity / delivery time upon request

Adaptiert auf .035 inch Führungsdraht
 adapté au fil-guide de .035 inch
 adapted to .035 inch guide wire

Verpackungseinheit: 1 Stück
 conditionnement : 1 unité
 box: 1 unit