


# sinus-Repo-Visual 6F

self-expanding nitinol stent system



the repositionable  
6F nitinol stent with  
closed-cell design for  
straight vessel segments



POWERFUL & REPOSITIONABLE

# sinus- Repo-Visual 6F





closed-cell design



repositionable stent



electro-polishing



braided sheath




radiopaque markers



atraumatic soft-tip



anti-jump technique



Bestell-Nr. – Länge des Einführbestecks référence – longueur système de pose order code – length of application device		Stent	Stent Länge / longueur / length
75 cm	120 cm	Ø / mm	mm
6306-6030	6306-7030	6	30
6306-6040	6306-7040	6	40
6306-6060	6306-7060	6	60
6306-6080	6306-7080*	6	80
6307-6020	6307-7020*	7	20
6307-6030	6307-7030*	7	30
6307-6040	6307-7040*	7	40
6307-6060	6307-7060*	7	60
6307-6080	6307-7080*	7	80
6308-6020	6308-7020*	8	20
6308-6030	6308-7030*	8	30
6308-6040	6308-7040*	8	40
6308-6060	6308-7060*	8	60
6308-6080	6308-7080*	8	80
6309-6030	6309-7030*	9	30
6309-6040	6309-7040*	9	40
6309-6060	6309-7060*	9	60
6309-6080	6309-7080*	9	80
6310-6030	6310-7030*	10	30
6310-6040	6310-7040*	10	40
6310-6060	6310-7060*	10	60
6310-6080	6310-7080*	10	80
6312-6030	6312-7030*	12	30
6312-6040	6312-7040*	12	40
6312-6060	6312-7060*	12	60
6314-6040	6314-7040*	14	40
6314-6060	6314-7060*	14	60
6314-6080	6314-7080*	14	80
6314-6100	6314-7100*	14	100

\* minimum order quantity / delivery time upon request  
 \* Mindestbestellmenge / Lieferzeit auf Anfrage  
 \* quantité minimum de commande / délai de livraison sur demande

Verpackungseinheit: 1 Stück  
 conditionnement : 1 unité  
 box: 1 unit

## FEATURES

- 6F application device adapted to 0.035 inch guide wire
- Repositionable self-expanding nitinol stent
- Maximum radial force for straight vessel segments and highly calcified stenosis

## BESONDERHEITEN

- 6F Applikationssystem adaptiert auf 0,035 inch Führungsdraht
- Repositionierbarer selbstexpandierbarer Nitinol Stent
- Höchste Radialkraft für gerade Gefäßabschnitte und stark kalkifizierte Stenosen

## CARACTERISTIQUES

- Accès 6F sur guide 0.035 inch
- Stent autoexpandible en Nitinol repositionnable
- Force radiale maximale pour les parties rectilignes des vaisseaux et les sténoses fortement calcifiées

+ the perfect supplement

### NyloTrack-35

semi-compliant PTA-balloon catheter

### MARS

non-compliant PTA-balloon catheter



radiopaque markers



re-wrapping with dog-bone



radiopaque tip



kink-resistant shaft